

Sterile suction cannula

Instructions for use

01 Intended use

A cylindrical suction cannula used for aspiration with a suction device. Connect the tubing to a suction pump.

02 Preparation before use

1. Check the expiry date of sterility.
2. Do not use if the expiry date has passed or the packaging is damaged.
3. Understand how to use the product before use.

03 Operation

1. Open the packaging in a sterile area.
2. Connect to the suction device and suction tubing for use.
3. Dispose of used products according to the prescribed procedures.
4. Do not reuse a product that has been used.

04 After-use care

After use, dispose of the product in accordance with applicable procedures. Do not reuse.

05 Precautions

1. Check the expiry date of sterility.
2. Do not use if the expiry date has passed or the packaging is damaged.
3. Understand how to use the product before use.
4. This is a sterile, single-use medical device. Do not disinfect again or reuse.

06 Storage

1. Store in a hygienic location at room temperature, 1-30°C.
2. After use, dispose of the product in accordance with applicable procedures. Do not reuse.
3. Store away from direct sunlight.

07 Sterilization

Ethylene oxide (EO) gas sterilization